DCM013-11
Ed. 01/2015

TSH ELISA

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Ed. 01/2015

DCM013-11

DiaMetra S.r.l. Headquarter: Via Calabria 15,
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

Manufactory: Via Pozzuolo 14
06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: info@diametra.com

DCM013-11
Ed. 01/2015

TSH ELISA

for routine analysis

Direct immunoenzymatic determination of TSH in human serum or plasma

IVD

LOT

See external label

Σ = 96 tests

REF DKO013

INTENDED USE

Diametra TSH ELISA is an immunoenzymatic colorimetric method for quantitative determination of thyroid-stimulating hormone (TSH, thyrotropin) concentration in human serum or plasma. TSH ELISA kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Thyroid-stimulating hormone (TSH or thyrotropin) is a glycopolypeptide hormone synthesized and secreted by the anterior pituitary gland which regulates the endocrine function of the thyroid gland.

TSH consists of two subunits, the *alpha* and the *beta* subunit. The α subunit is identical to that of human chorionic gonadotropin (HCG), luteinising hormone (LH), follicle-stimulating hormone (FSH). The (*beta*) subunit is unique to TSH, and therefore determines its function.

TSH production is controlled by a Thyrotropin Releasing Hormone, (TRH), which is manufactured in the hypothalamus and transported to the pituitary gland, where it increases TSH production and release. Somatostatin is also produced by the hypothalamus, and has an opposite effect on the pituitary production of TSH, decreasing or inhibiting its release.

TSH stimulates the thyroid gland to secrete the hormones thyroxine (T4) and triiodothyronine (T3).

Release of TSH is regulated by the circulating free fraction of thyroid hormones in the blood. TSH levels are depressed when peripheral concentrations of the free fraction of thyroid hormones are high. Conversely, TSH levels are high when peripheral concentrations of thyroid hormones are low. This effect creates a regulatory negative feedback loop.

TSH levels are tested in the blood of patients suspected of suffering from excess (hyperthyroidism), or deficiency (hypothyroidism) of thyroid hormone. Generally, a normal range for TSH is between 0.3 and 3.0 mIU/mL, but the interpretation depends also on what the blood levels of thyroid hormones (T3 and T4) are.

Higher than normal levels of TSH combined with high levels of thyroid hormone (T3 and T4) may indicate dysfunction of the hypothalamus and pituitary gland. In these case, a high TSH is often produced by a benign

tumour of the pituitary (adenoma). Conversely, low levels of TSH, while blood levels of T3 and T4 are also low, indicates abnormally low function of the pituitary, known as hypopituitarism.

On the other hand, due to the negative feedback described above, abnormally high levels of Thyroid hormone, due to overproduction in the thyroid, results in low TSH levels. This occurs in diseases such as hyperthyroidism or Grave's disease. Conversely, an underproduction of T3 and T4 caused by diseases such as congenital hypothyroidism (cretinism), hypothyroidism or thyroid hormone resistance, gives rise to an increase in the measured TSH.

Clearly both TSH and T3 and T4 should be measured to ascertain where a specific thyroid disfunction is caused by primary pituitary or by a primary thyroid disease.


2. PRINCIPLE

The essential reagents required for an immunoenzymatic assay include high affinity and specificity antibodies (enzyme and immobilised) with different and distinct epitope recognition, in excess, and native antigen.

In this method, TSH calibrators, patient specimens and/or controls containing the native antigen are first added to streptavidin coated wells. Biotinilated monoclonal and enzyme labeled antibodies are added and the reactants mixed : these antibodies have high affinity and specificity and are detected against distinct and different epitopes of TSH.

Reaction between the various TSH antibodies and native TSH occurs in the microwells without competition or steric hindrance forming a soluble sandwich complex.

The interaction is illustrated by the following equation:

BtAb(m) = Biotinilated Monoclonal Antibody (Excess Quantity)

AgTSH = Native Antigen (Variable Quantity)

EAb(p) = Enzyme labeled Antibody (Excess Quantity)

EAb(p)-AgTSH-BtAb(m) = Antigen-Antibodies Sandwich Complex

Ka = Rate Constant of Association

K-a = Rate Constant of Dissociation

Simultaneously, the complex is deposited to the well though the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:

EAb(p)-AgTSH-BtAb(m) + Streptavidina CW. \Rightarrow Immobilized Complex

Streptavidin CW. = Streptavidin immobilized on well
Immobilized Complex = Antibodies-Antigen sandwich bound.

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody bound fraction is directly proportional to the native antigen concentration. The activity of the enzyme present on the surface of the well is quantified by reaction with a suitable substrate to produce colour. By utilizing several different calibrators of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. Calibrators (7 vials, 1 mL each)

CAL0	REF DCE002/1306-0
CAL1	REF DCE002/1307-0
CAL2	REF DCE002/1308-0
CAL3	REF DCE002/1309-0
CAL4	REF DCE002/1310-0
CAL5	REF DCE002/1311-0
CAL6	REF DCE002/1312-0

2. Control (1 vial, 1 mL)

Concentration of Control is indicated on the Certificate of Analysis

REF DCE045/1303-0

3. Conjugate (1 vial, 13 mL)

Goat antibody anti TSH conjugated with horseradish peroxidase HRP; mouse antibody anti TSH biotinylated

REF DCE002/1302-0

4. Coated Microplate (1 breakable microplate)

Microplate coated with streptavidin

REF DCE002/1303-0

5. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB 0.26 g/L (avoid any skin contact)

REF DCE004-0

6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact)

REF DCE005-0

7. 50X Conc. Wash Solution (1 vial, 20 mL)

NaCl 45 g/L, Tween-20 55 g/L

REF DCE006-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm).

Note

Store all reagents between 20-8°C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close immediately after use. Once opened, the microplate is stable until expiry date of the kit.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300^R as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of TSH from 0.2 mIU/L to 20.0 mIU/L.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are

needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate

- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrators (C₀...C₆)

The Calibrators are ready to use, are calibrated against WHO 2nd IRP 80/558 and have the following concentrations of TSH:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆
mIU/L	0	0.2	0.5	2.5	5.0	10	20

The Calibrators are stable until the expiry date printed on the label. Once opened, the calibrators are stable six months at 2-8°C.

6.2. Preparation of Wash Solution

Dilute the contents of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

6.3. Preparation of the Sample

TSH assay can be performed in human serum or plasma.

Specimens may be stored at 2-8°C for maximum two days. For longer storage, specimens should be frozen. Frozen specimens should be well mixed after thawing, and before assay. Avoid repeated freezing and thawing.

Serum samples with TSH concentrations greater than 20mIU/L should be diluted 1:2 with C₀. Correct the result using an appropriate dilution factor.

The Control is ready to use.

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.

- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₆), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C ₀ -C ₆	50 µL		
Sample/Control		50 µL	
Conjugate	100 µL	100 µL	
Cover the plate and incubate for 60 minutes at room temperature (22±28°C). To increase the sensibility read the absorbance after 120 min of incubation. Remove the content from each well; wash the wells 6 times with 300 µL of diluted Wash Solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.			
TMB Substrate	100 µL	100 µL	100 µL
Cover the plate and incubate for 20 minutes at room temperature (22±28°C) in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the plate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of TSH for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

8.1. Mean Absorbance

Calculate the mean of the absorbance (Em) for each point of the calibration curve (C₀-C₆) and of each sample.

8.2. Calibration curve

Plot the mean value of absorbance (Em) of the calibrators (C₀-C₆) against concentration. Draw the best-fit curve through the plotted points. (es:Cubic Spline or Four Parameter Logistic).

8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in mIU/L.

9. REFERENCE VALUES

Serum samples of apparently healthy women and men were assayed using the TSH ELISA test, with the following results:

	Mean (mIU/L)	Range (mIU/L)
TSH	1.85	0.39 – 6.16

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

10.1.1. Intra Assay Variation

Within run variation was determined by replicate measurements (16x) of two different control sera in one assay. The within assay variability is 4.6%.

10.1.2. Inter Assay Variation

Between run variation was determined by replicate measurements (12x) of three different control sera in different lots of kit. The between assay variability is 10.8%.

10.2. Accuracy

The recovery of 0.25 – 2.5 – 10 mIU/L of TSH added to sample gave an average value (±SD) of 85.85% ± 2.62% with reference to the original concentrations.

10.3. Sensitivity

The lowest detectable concentration of TSH that can be distinguished from the Calibrator 0 is 0.01 mIU/L with a confidence limit of 95%.

10.4. Specificity

The cross-reactivity of the thyrotropin ELISA method to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated

by the ratio between the dose of the interfering substance and the amount of TSH needed to obtain the same absorbance.

Thyrotropin (hTSH)	100 %
Thyrotropin - alpha (hTSH-)	9.3 %
Thyrotropin - beta (hTSH-)	152 %
Follicle-Stimulating Hormone (hFSH)	1.0 %
Luteinising hormone (hLH)	1.0 %
Chorionic Gonadotropin (hCG)	<0.1 %

10.5. Correlation with CLIA

Diametra TSH ELISA was compared with two others commercially available TSH assays. 54 serum samples were tested.

The two linear regression curves were calculated:

$$(TSH \text{ Diametra}) = 0.96 * (TSH \text{ CLIA}) + 0.06$$

$$r^2 = 0.972$$

$$(TSH \text{ Diametra}) = 0.87 * (TSH \text{ Elisa}) + 0.29$$

$$r^2 = 0.971$$

10.6. Hook Effect

Diametra TSH ELISA, a competitive enzyme immunoassay, shows no Hook Effect up to 500 mIU/L.

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

- Hopton M. R, et al *Clinical Chemistry*, 32, 691 (1986)
- Caldwell, G. et al. *Lancet*, I, 1117, (1985)
- Young, D. S, et al *Clinical Chemistry* 21, 3660 (1975)
- Spencer, CA, et al *Clinical Chemistry* 41, 367 (1995)
- Beck-Peccoz P, et al *Eur. J. Endocrinol* 131: 331-340 (1994)
- Bravemann, L. E, *Clinical Chemistry* 42: 174-181 (1996)
- Fisher, D. A., *Clinical Chemistry* 42: 135-139 (1996)

Ed. 01/2015

DCM013-11

DiaMetra S.r.l. Headquarter: Via Calabria 15 –
20090 SEGRATE (MI) Italy
Tel. +39.02.2139184;
Fax +39.02.2133354.

Manufactory: Via Pozzuolo 14, 06038 SPELLO (PG)
Italy
Tel. +39.0742.24851;
Fax +39.0742.316197
E-mail: info@diametra.com

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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