

订购热线: 4008-898-798

# Anti-SSX1 antibody

**Cat. No.** ml123233

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-SSX1 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Fusion protein of human SSX1

ReactivityHumanContent0.4 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol SSX1

Full name synovial sarcoma, X breakpoint 1

Synonyms SSRC; CT5.1 Swissprot Q16384

### **Target Background**

The product of this gene belongs to the family of highly homologous synovial sarcoma X (SSX) breakpoint proteins. These proteins may function as transcriptional repressors. They are also capable of eliciting spontaneous humoral and cellular immune responses in cancer patients, and are potentially useful targets in cancer vaccine-based immunotherapy.



订购热线: 4008-898-798

## **Applications**

## Western blotting

Predicted band size:22 kDa Positive control:K562 cells

Recommended dilution: 500-2000

Gel: 12%SDS-PAGE


Lysate: 40 µg Lane: K562 cells

Primary antibody: ml123233(SSX1 Antibody) at dilution

1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



### **ELISA**

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn