

# Anti-KDM4C antibody

 Cat. No.
 ml161810

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-KDM4C rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human KDM4C
Reactivity	Human, Mouse
Content	1 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	KDM4C
Full name	lysine (K)-specific demethylase 4C

GASC1; JHDM3C; JMJD2C; TDRD14C; bA146B14.1 Q9H3R0

#### **Target Background**

Synonyms

Swissprot


This gene is a member of the Jumonji domain 2 (JMJD2) family and encodes a protein with one JmjC domain, one JmjN domain, two PHD-type zinc fingers, and two Tudor domains. This nuclear protein functions as a trimethylation-specific demethylase, converting specific trimethylated histone residues to the dimethylated form. Chromosomal aberrations and increased transcriptional expression of this gene are associated with esophageal squamous cell carcinoma. Alternative splicing results in multiple transcript variants.



## 订购热线: 4008-898-798

Applications Western blotting Predicted band size:120 kDa Positive control: Mouse heart and brain tissue Recommended dilution: 500-2000

Gel: 6%SDS-PAGE Lysate: 40 µg Lane 1-2: Mouse heart tissue, Mouse brain tissue Primary antibody: ml161810(KDM4C Antibody) at dilution 1/800 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution Exposure time: 1 minute



#### ELISA

Recommended dilution:2000-5000

Wit produce 联系电话: 4008-898-798, 021-61725725

### 联系QQ: 2881505695, 2881505696

- 邮箱: mlbio\_cn@yeah.net
- 网址: www.mlbio.cn