

Anti-HTATIP2 antibody

Cat. No.	ml123355
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-HTATIP2 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human HTATIP2
Reactivity	Human, Mouse
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	HTATIP2
Full name	HIV-1 Tat interactive protein 2, 30kDa
Synonyms	CC3; TIP30; SDR44U1
Swissprot	Q9BUP3

Target Background

TIP30 (30 kDa HIV-1 TAT-interacting protein), also known as HTATIP2 (HIV-1 TAT-interactive protein 2) or CC3, is a ubiquitously expressed protein that functions as an oxidoreductase. It acts as a tumor suppressor participating in metabolic suppression, inhibition of angiogenesis and induction of apoptosis.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:27 kDa

Positive control:A431 cells

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: A431 cells

Primary antibody: ml123355(HTATIP2 Antibody) at dilution
1/800

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds



ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn