

Anti-XRN2 antibody

Cat. No.	ml124080
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-XRN2 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human XRN2
Reactivity	Human, Mouse
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	XRN2
Full name	5'-3' exoribonuclease 2
Synonyms	
Swissprot	Q9H0D6

Target Background

This gene shares similarity with the mouse Dhml and the yeast dhp1 gene. The yeast gene is involved in homologous recombination and RNA metabolism, such as RNA synthesis and RNA trafficking. Complementation studies show that Dhml has a similar function in mouse as dhp1. The function of the human gene has not yet been determined. Transcript variants encoding different isoforms have been noted for this gene; however, their full-length nature is not known.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:109 kDa

Positive control:Hela, NIH/3T3, A172 and HEPG2 cell lysates

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-4: Hela, NIH/3T3, A172 and HEPG2 cell lysates

Primary antibody: ml124080(XRN2 Antibody) at dilution 1/350

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn