

Anti-NADK2 antibody

Cat. No.	ml124363
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-NADK2 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human NADK2
Reactivity	Human, Rat
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	NADK2
Full name	NAD kinase 2, mitochondrial
Synonyms	DECRD; MNADK; NADKD1; C5orf33
Swissprot	Q4G0N4

Target Background

This gene encodes a mitochondrial kinase that catalyzes the phosphorylation of NAD to yield NADP. Mutations in this gene result in 2,4-dienoyl-CoA reductase deficiency. Alternative splicing results in multiple transcript variants.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:49 kDa

Positive control:HEPG2, HeLa and A549 cell lysates

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 μ g

Lane 1-3: HEPG2, HeLa and A549 cell lysates

Primary antibody: ml124363(NADK2 Antibody) at dilution 1/600

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 60 seconds



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn