

Anti-TNNI3 antibody

Cat. No.	ml220341
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-TNNI3 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human TNNI3
Reactivity	Human, Mouse, Rat
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	TNNI3
Full name	troponin I type 3 (cardiac)
Synonyms	CMH7; RCM1; cTnI; CMD2A; TNNC1; CMD1FF
Swissprot	P19429

Target Background

Troponin I (TnI), along with troponin T (TnT) and troponin C (TnC), is one of 3 subunits that form the troponin complex of the thin filaments of striated muscle. TnI is the inhibitory subunit; blocking actin-myosin interactions and thereby mediating striated muscle relaxation. The TnI subfamily contains three genes: TnI-skeletal-fast-twitch, TnI-skeletal-slow-twitch, and TnI-cardiac. This gene encodes the TnI-cardiac protein and is exclusively expressed in cardiac muscle tissues. Mutations in this gene cause familial hypertrophic cardiomyopathy type 7 (CMH7) and familial restrictive cardiomyopathy (RCM).

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 24 kDa

Positive control: Mouse heart tissue

Recommended dilution: 200-1000

Gel: 12% SDS-PAGE


Lysate: 40 µg

Lane: Mouse heart tissue

Primary antibody: ml220341 (TNNI3 Antibody) at dilution 1/175

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn