

Anti-MCM9 antibody

Cat. No.	ml220964
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-MCM9 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human MCM9
Reactivity	Human, Mouse, Rat
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	MCM9
Full name	minichromosome maintenance 9 homologous recombination repair factor
Synonyms	ODG4; MCMDC1; C6orf61; dJ329L24.1; dJ329L24.3
Swissprot	Q9NXL9

Target Background

The protein encoded by this gene is a member of the mini-chromosome maintenance (MCM) protein family that are essential for the initiation of eukaryotic genome replication. Binding of this protein to chromatin has been shown to be a prerequisite for recruiting the MCM2-7 helicase to DNA replication origins. This protein also binds, and is a positive regulator of, the chromatin licensing and DNA replication factor 1, CDT1.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:127 kDa

Positive control:K562 and 231 cell lysates

Recommended dilution: 500-2000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-2: K562 and 231 cell lysates

Primary antibody: ml220964(MCM9 Antibody) at dilution 1/800

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn