

# Anti-CSNK2A1 antibody

**Cat. No.** ml220424

Package 25  $\mu$ l/100  $\mu$ l/200  $\mu$ l

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-CSNK2A1 rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Fusion protein of human CSNK2A1

Reactivity Human, Mouse, Rat

Content 0.3 mg/ml

Host species Rabbit

Ig class Immunogen-specific rabbit IgG

Purification Antigen affinity purification

**Target information** 

Symbol CSNK2A1

Full name casein kinase 2, alpha 1 polypeptide



Synonyms CKII; CK2A1; CSNK2A3

Swissprot P68400

#### **Target Background**

Casein kinase II is a serine/threonine protein kinase that phosphorylates acidic proteins such as casein. It is involved in various cellular processes, including cell cycle control, apoptosis, and circadian rhythms. The kinase exists as a tetramer and is composed of an alpha, an alpha-prime, and two beta subunits. The alpha subunits contain the catalytic activity while the beta subunits undergo autophosphorylation. The protein encoded by this gene represents the alpha subunit. While this gene is found on chromosome 20, a related transcribed pseudogene is found on chromosome 11. Three transcript variants encoding two different proteins have been found for this gene.




## **Applications**


## Immunohistochemistry

Predicted cell location: Cytoplasm and Nucleus

Positive control: Human liver cancer

Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml220424(CSNK2A1 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

#### Western blotting

Predicted band size:45 kDa

Positive control:K562 and 293T cells, Jurkat cells and mouse brain tissue

Good elisakii produceri

Recommended dilution: 500-2000



Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-4: K562 cells, 293T cells, Jurkat cells, mouse brain tissue

Primary antibody: ml220424(CSNK2A1 Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



## **ELISA**

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn