

# 订购热线: 4008-898-798

## Anti-CNGA3 antibody

Cat. No. ml220446 Package 25 µl/100 µl/200 µl Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-CNGA3 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human CNGA3
Reactivity	Human, Mouse
Content	0.7 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	CNGA3
Full name	cyclic nucleotide gated channel alpha 3

30000

Symbol	CNGA3
Full name	cyclic nucleotide gated channel alpha 3
Synonyms	CNG3; ACHM2; CCNC1; CCNCa; CNCG3; CCNCalpha
Swissprot	Q16281
Target Background	

#### **Target Background**


This gene encodes a member of the cyclic nucleotide-gated cation channel protein family which is required for normal vision and olfactory signal transduction. Mutations in this gene are associated with achromatopsia (rod monochromacy) and color blindness. Two alternatively spliced transcripts encoding different isoforms have been described.



### 订购热线: 4008-898-798

Applications Western blotting Predicted band size:79 kDa Positive control: Mouse kidney and brain tissue Recommended dilution: 500-2000

Gel: 6%SDS-PAGE Lysate: 40 µg Lane 1-2: Mouse kidney tissue, Mouse brain tissue Primary antibody: ml220446(CNGA3 Antibody) at dilution 1/500 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution Exposure time: 1 minute



#### ELISA

Recommended dilution:2000-5000

Wit produce 联系电话: 4008-898-798, 021-61725725

### 联系QQ: 2881505695, 2881505696

- 邮箱: mlbio\_cn@yeah.net
- 网址: www.mlbio.cn