

## Anti-TNFRSF12A antibody

<b>Cat. No.</b>	ml220536
<b>Package</b>	25 µl/100 µl/200 µl
<b>Storage</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol

### Product overview

<b>Description</b>	Anti-TNFRSF12A rabbit polyclonal antibody
<b>Applications</b>	ELISA, WB
<b>Immunogen</b>	Fusion protein of human TNFRSF12A
<b>Reactivity</b>	Human, Mouse
<b>Content</b>	0.1 mg/ml
<b>Host species</b>	Rabbit
<b>Ig class</b>	Immunogen-specific rabbit IgG
<b>Purification</b>	Antigen affinity purification

### Target information

<b>Symbol</b>	TNFRSF12A
<b>Full name</b>	Tumor necrosis factor receptor superfamily, member 12A
<b>Synonyms</b>	FN14; CD266; TWEAKR
<b>Swissprot</b>	Q9NP84

### Target Background

Tumor necrosis factor receptor superfamily member 12A also known as the Tweak-receptor (TweakR) is a protein that in humans is encoded by the TNFRSF12A gene. Receptor for TNFSF12/TWEAK. Weak inducer of apoptosis in some cell types. Promotes angiogenesis and the proliferation of endothelial cells. May modulate cellular adhesion to matrix proteins. Highly expressed in heart, placenta and kidney. Intermediate expression in lung, skeletal muscle and pancreas.

订购热线: 4008-898-798

**Applications**

**Western blotting**

Predicted band size: 14 kDa

Positive control: 293T cells

Recommended dilution: 500-2000

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane: 293T cells

Primary antibody: ml220536(TNFRSF12A Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



**ELISA**

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn