

Anti-TIMP4 antibody

Cat. No.	ml221340
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-TIMP4 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human TIMP4
Reactivity	Human, Mouse, Rat
Content	0.1 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	TIMP4
Full name	TIMP metalloproteinase inhibitor 4
Synonyms	
Swissprot	Q99727

Target Background

This gene belongs to the TIMP gene family. The proteins encoded by this gene family are inhibitors of the matrix metalloproteinases, a group of peptidases involved in degradation of the extracellular matrix. The secreted, netrin domain-containing protein encoded by this gene is involved in regulation of platelet aggregation and recruitment and may play role in hormonal regulation and endometrial tissue remodeling.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 26 kDa

Positive control: Human brain tissue, PC3 and hela cells

Recommended dilution: 200-1000

Gel: 8+12% SDS-PAGE


Lysate: 50 µg

Lane 1-3: Human brain tissue, PC3 cells, hela cells

Primary antibody: ml221340 (TIMP4 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn