

Anti-IL18RAP antibody

Cat. No.	ml221404
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-IL18RAP rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human IL18RAP
Reactivity	Human, Mouse
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	IL18RAP
Full name	Interleukin 18 receptor accessory protein
Synonyms	ACPL; CD218b; IL18RB; CDw218b
Swissprot	O95256

Target Background

Interleukin 18 receptor accessory protein, also known as IL18RAP and CDw218b (cluster of differentiationw218b), is a human gene. The protein encoded by this gene is an accessory subunit of the heterodimeric receptor for IL18. This protein enhances the IL18 binding activity of IL18R1 (IL1RRP), a ligand binding subunit of IL18 receptor. The coexpression of IL18R1 and this protein is required for the activation of NF-κB and MAPK8 (JNK) in response to IL18.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 68 kDa

Positive control: RAW264.7 cells and human liver cancer tissue, SP20 cells

Recommended dilution: 500-2000

Gel: 10% SDS-PAGE


Lysate: 40 µg

Lane 1-3: RAW264.7 cells, human liver cancer tissue, SP20 cell

Primary antibody: ml221404(IL18RAP Antibody) at dilution 1/650

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn