

兔抗 CITED1 多克隆抗体

中文名称：兔抗 CITED1 多克隆抗体

英文名称：Anti-CITED1 rabbit polyclonal antibody

别名：Cbp/p300 interacting transactivator with Glu/Asp rich carboxy-terminal domain 1; MSG1

相关类别：一抗

抗原：CITED1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the CREB-binding protein/p300-interacting transactivator with Asp/Glu-rich C-terminal domain (CITED) family of proteins. The encoded protein, also known as melanocyte-specific gene 1, may function as a transcriptional coactivator and may play a role in pigmentation of melanocytes. Alternatively spliced transcript variants have been described.
Applications:	ELISA, WB, IHC
Name of antibody:	CITED1
Immunogen:	Synthetic peptide of human CITED1

Full name:	Cbp/p300 interacting transactivator with Glu/Asp rich carboxy-terminal domain 1
Synonyms:	MSG1
SwissProt:	Q99966
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	20 kDa
WB Positive control:	Human testis tissue lysate
WB Recommended dilution:	200-1000



