

## 兔抗 CKAP2 多克隆抗体

中文名称：兔抗 CKAP2 多克隆抗体

英文名称： Anti-CKAP2 rabbit polyclonal antibody

别名： LB1; TMAP; se20-10

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： CKAP2

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a cytoskeleton-associated protein that stabilizes microtubules and plays a role in the regulation of cell division. The encoded protein is itself regulated through phosphorylation at multiple serine and threonine residues. There is a pseudogene of this gene on chromosome 14. Alternative splicing results in multiple transcript variations.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CKAP2
<b>Immunogen:</b>	Fusion protein of human CKAP2
<b>Full name:</b>	cytoskeleton associated protein 2

<b>Synonyms :</b>	LB1; TMAP; se20-10
<b>SwissProt:</b>	Q8WWK9
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human colon cancer and human tonsil
<b>IHC Recommend dilution:</b>	50-200

