

兔抗 ANKRA2 多克隆抗体

中文名称：兔抗 ANKRA2 多克隆抗体

英文名称： Anti-ANKRA2 rabbit polyclonal antibody

别名： ANKRA

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ANKRA2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This protein may facilitate endocytosis by linking megalin to components of the cytoskeleton or endocytic machinery. ANKRA2 is a megalin- and BKCa potassium channel-interacting factor. It is a protein known to interact with HDAC4 and HDAC5. ANKRA2 is important for transcriptional repression by AhRR.
Applications:	ELISA, IHC
Name of antibody:	ANKRA2
Immunogen:	Synthetic peptide of human ANKRA2
Full name:	ankyrin repeat, family A (RFXANK-like), 2
Synonyms :	ANKRA

SwissProt:	Q9H9E1
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer and Human lung cancer
IHC Recommend dilution:	25-100

