

兔抗 CLASP1 多克隆抗体

中文名称：兔抗 CLASP1 多克隆抗体

英文名称：Anti-CLASP1 rabbit polyclonal antibody

别名：cytoplasmic linker associated protein 1; MAST1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CLASP1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	CLASPs, such as CLASP1, are nonmotor microtubule-associated proteins that interact with CLIPs (e.g., CLIP170; MI M 179838). CLASP1 is involved in the regulation of microtubule dynamics at the kinetochore and throughout the spindle (Maiato et al., 2003 [PubMed 12837247]).
Applications:	ELISA, IHC
Name of antibody:	CLASP1
Immunogen:	Fusion protein of human CLASP1
Full name:	cytoplasmic linker associated protein 1
Synonyms:	MAST1
SwissProt:	Q7Z460

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	20-100

