

兔抗 CLCNKA 多克隆抗体

中文名称：兔抗 CLCNKA 多克隆抗体

英文名称： Anti-CLCNKA rabbit polyclonal antibody

别名： chloride voltage-gated channel Ka; CLCK1; ClC-K1; hClC-Ka

相关类别： 一抗

抗原： CLCNKA

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a member of the CLC family of voltage-gated chloride channels. The encoded protein is predicted to have 12 transmembrane domains, and requires a beta subunit called barttin to form a functional channel. It is thought to function in salt reabsorption in the kidney and potassium recycling in the inner ear. The gene is highly similar to CLCNKB, which is located 10 kb downstream from this gene. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	CLCNKA

Immunogen:	Synthetic peptide of human CLCNKA
Full name:	chloride voltage-gated channel Ka
Synonyms:	CLCK1; CIC-K1; hCIC-Ka
SwissProt:	P51800
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	20-100

