

兔抗 CLDN3 多克隆抗体

中文名称: 兔抗 CLDN3 多克隆抗体

英文名称: Anti-CLDN3 rabbit polyclonal antibody

别 名: RVP1, HRVP1, C7orf1, CPE-R2, CPETR3

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: CLDN3

反应种属: Human

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

IHC Recommend dilution:	25-100
IHC positive control:	Human breast cancer and Human colon cancer
ELISA Recommended dilution:	1000-10000
SwissProt:	O15551
Synonyms:	RVP1, HRVP1, C7orf1, CPE-R2, CPETR3
Full name:	claudin 3
Immunogen:	Synthetic peptide of human CLDN3
Name of antibody:	CLDN3
Applications:	ELISA, IHC
Background:	Tight junctions represent one mode of cell-to-cell adhesi on in epithelial or endothelial cell sheets, forming contin



uous seals around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. These junctions are comprised of s ets of continuous networking strands in the outwardly fa cing cytoplasmic leaflet, with complementary grooves in t he inwardly facing extracytoplasmic leaflet. The protein e ncoded by this intronless gene, a member of the claudin family, is an integral membrane protein and a componen t of tight junction strands. It is also a low-affinity recepto r for Clostridium perfringens enterotoxin, and shares aa s equence similarity with a putative apoptosis-related protein found in rat.

