

兔抗 CLDN5 多克隆抗体

中文名称：兔抗 CLDN5 多克隆抗体

英文名称： Anti-CLDN5 rabbit polyclonal antibody

别名： AWAL, BEC1, TMVCF, CPETRL1

相关类别： 一抗

抗原： CLDN5

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets. Mutations in this gene have been found in patients with velocardiofacial syndrome. Alternatively spliced transcript variants encoding the same protein have been found for this gene.
Applications:	ELISA, IHC

Name of antibody:	CLDN5
Immunogen:	Synthetic peptide of human CLDN5
Full name:	claudin 5
Synonyms :	AWAL, BEC1, TMVCF, CPETRL1
SwissProt:	O00501
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	5-20

