

兔抗 ANXA1(Ab-21) 多克隆抗体

中文名称：兔抗 ANXA1(Ab-21) 多克隆抗体

英文名称： Anti-ANXA1(Ab-21) rabbit polyclonal antibody

别名： ANX1; LPC1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： ANXA1(Ab-21)

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis. This protein regulates phospholipase A2 activity. It seems to bind from two to four calcium ions with high affinity.
Applications:	WB
Name of antibody:	ANXA1(Ab-21)
Immunogen:	Synthesized non-phosphopeptide derived from human ANXA1 around the phosphorylation site of tyrosine 21 (Q-E-Y(p)-V-Q).
Full name:	annexin A1
Synonyms :	ANX1; LPC1

SwissProt:	P04083
WB Predicted band size:	39 kDa
WB Positive control:	293 cells lysate
WB Recommended dilution:	500-3000

