

兔抗 CLDN6 多克隆抗体

中文名称: 兔抗 CLDN6 多克隆抗体

英文名称: Anti-CLDN6 rabbit polyclonal antibody

抗 原: CLDN6

储 存: 冷冻 (-20℃)

宿 主: Rabbit

反应种属: Human, Mouse

标 记 物: Unconjugate

相关类别: 一抗

克隆类型: rabbit polyclonal

技术规格

Background:

Tight junctions represent one mode of cell-to-cell adhesion in epithelial or endothelial cell sheets, forming continuous seals around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. These junctions are comprised of sets of continuous networking strands in the outwardly facing cytoplasmic leaflet, with complementary grooves in the inwardly facing extracytoplasmic leaflet. This gene encodes a component of tight junction strands, which is a member of the claudin family. The protein is an integral membrane protein and is one of the entry cofactors for hepatitis C virus. The gene met

	ylation may be involved in esophageal tumorigenesis. This gene is adjacent to another family member CLDN 9 on chromosome 16.
Applications:	ELISA, IHC
Name of antibody:	CLDN6
Immunogen:	Synthetic peptide of human CLDN6
Full name:	claudin 6
SwissProt:	P56747
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human esophagus cancer
IHC Recommend dilution:	25-100



