

## 兔抗 CLDN7(Ab-210) 多克隆抗体

中文名称：兔抗 CLDN7(Ab-210) 多克隆抗体

英文名称：Anti-CLDN7(Ab-210) rabbit polyclonal antibody

别名：CLDN-7; CEPTRL2; CPETRL2; Hs.84359; claudin-

相关类别：一抗

抗原：CLDN7(Ab-210)

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Plays a major role in tight junction-specific obliteration of the intercellular space By similarity.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	CLDN7(Ab-210)
<b>Immunogen:</b>	Synthesized non-phosphopeptide derived from human Claudin 7 around the phosphorylation site of tyrosine 210 (S-K-E-Y(p)-V).
<b>Full name:</b>	claudin 7
<b>Synonyms :</b>	CLDN-7; CEPTRL2; CPETRL2; Hs.84359; claudin-
<b>SwissProt:</b>	O95471
<b>IHC positive control:</b>	Human prostate carcinoma tissue

<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	22 kDa
<b>WB Positive control:</b>	Rat liver tissue lysate
<b>WB Recommended dilution:</b>	500-3000



