

兔抗 CLDN8 多克隆抗体

中文名称: 兔抗 CLDN8 多克隆抗体

英文名称: Anti-CLDN8 rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: CLDN8

反应种属: Human

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

This gene encodes a member of the claudin family. Cla udins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epit helial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. This protein plays important roles in the paracellular cation barrier of the distal renal tubule, and in the paracellular barrier to prevent sodium back-leakage in distal colon. Differential expression of this gene has been observed in colorectal carcinoma and renal cell tumors, and along with claudin-7, is an immunohistochemical mar



	ker for the differential diagnosis of chromophobe renal cell carcinoma and renal oncocytoma.
Applications:	ELISA, IHC
Name of antibody:	CLDN8
Immunogen:	Synthetic peptide of human CLDN8
Full name:	claudin 8
SwissProt:	P56748
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human gastric cancer
IHC Recommend dilution:	50-200

