

## 兔抗 ANXA2 多克隆抗体

- 中文名称： 兔抗 ANXA2 多克隆抗体
- 英文名称： Anti-ANXA2 rabbit polyclonal antibody
- 别名： P36, ANX2, LIP2, LPC2, CAL1H, LPC2D, ANX2L4, PAP-IV
- 相关类别： 一抗
- 抗原： ANXA2
- 储存： 冷冻（-20℃）
- 宿主： Rabbit
- 反应种属： Human, Mouse, Rat
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. This protein functions as an autocrine factor which heightens osteoclast formation and bone resorption. This gene has three pseudogenes located on chromosomes 4, 9 and 10, respectively. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, WB, IHC

<b>Name of antibody:</b>	ANXA2
<b>Immunogen:</b>	Fusion protein of human ANXA2
<b>Full name:</b>	Annexin A2
<b>Synonyms :</b>	P36, ANX2, LIP2, LPC2, CAL1H, LPC2D, ANX2L4, PAP-IV
<b>SwissProt:</b>	P07355
<b>ELISA Recommended dilution:</b>	1000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	38 kDa
<b>WB Positive control:</b>	HepG2 cells
<b>WB Recommended dilution:</b>	1000-5000



