

兔抗 ANXA8 多克隆抗体

中文名称：兔抗 ANXA8 多克隆抗体

英文名称： Anti-ANXA8 rabbit polyclonal antibody

别名： ANX8

抗原： ANXA8

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the annexin family of evolutionarily conserved Ca ²⁺ and phospholipid binding proteins. The encoded protein may function as an anticoagulant that indirectly inhibits the thromboplastin-specific complex. Overexpression of this gene has been associated with acute myelocytic leukemia. A highly similar duplicated copy of this gene is found in close proximity on the long arm of chromosome 10.
Applications:	ELISA, IHC
Name of antibody:	ANXA8
Immunogen:	Fusion protein of human ANXA8

Full name:	annexin A8
Synonyms :	ANX8
SwissProt:	P13928
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human brain and human lung cancer
IHC Recommend dilution:	25-100

