

兔抗 ANXA10 多克隆抗体

中文名称：兔抗 ANXA10 多克隆抗体

英文名称：Anti-ANXA10 rabbit polyclonal antibody

别名：ANX14

相关类别：一抗

储存：冷冻（-20℃）

抗原：ANXA10

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cell growth and in signal transduction pathways. The function of this gene has not yet been determined.
Applications:	ELISA, WB, IHC
Name of antibody:	ANXA10
Immunogen:	Full length fusion protein
Full name:	annexin A10
Synonyms :	ANX14
SwissProt:	Q9UJ72

ELISA Recommended dilution:	1000-2000
IHC positive control:	Human breast cancer and Human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	37 kDa
WB Positive control:	A431 cell
WB Recommended dilution:	200-1000



