

兔抗 CLDN15 多克隆抗体

中文名称: 兔抗 CLDN15 多克隆抗体

英文名称: Anti-CLDN15 rabbit polyclonal antibody

抗 原: CLDN15

储 存:冷冻(-20℃)

宿 主: Rabbit

反应种属: Human

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the claudin family. Clau dins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. Alternatively spliced transcript variants encoding the same protein have been found for this gene.
Applications:	ELISA, WB
Name of antibody:	CLDN15
Immunogen:	Synthetic peptide of human CLDN15
Full name:	claudin 15



SwissProt:	P56746
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	24 kDa
WB Positive control:	K562, 231 and PC3 cell
WB Recommended dilution:	200-1000

