

兔抗 ANXA13 多克隆抗体

中文名称：兔抗 ANXA13 多克隆抗体

英文名称： Anti-ANXA13 rabbit polyclonal antibody

别名： ISA; ANX13

抗原： ANXA13

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The specific function of this gene has not yet been determined; however, it is associated with the plasma membrane of undifferentiated, proliferating endothelial cells and differentiated villus enterocytes. Alternatively spliced transcript variants encoding different isoforms have been identified.
Applications:	ELISA, WB, IHC

Name of antibody:	ANXA13
Immunogen:	Fusion protein of human ANXA13
Full name:	annexin A13
Synonyms:	ISA; ANX13
SwissProt:	P27216
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human prostate cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	35 kDa
WB Positive control:	Mouse small intestines tissue lysate
WB Recommended dilution:	500-2000



