

## 兔抗 AOA 多克隆抗体

中文名称：兔抗 AOA 多克隆抗体

英文名称： Anti-AOA rabbit polyclonal antibody

别名： acyloxyacyl hydrolase

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： AOA

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This locus encodes both the light and heavy subunits of a cyloxyacyl hydrolase. The encoded enzyme catalyzes the hydrolysis of acyloxylacyl-linked fatty acyl chains from bacterial lipopolysaccharides, effectively detoxifying these molecules. The encoded protein may play a role in modulating host inflammatory response to gram-negative bacteria. Alternatively spliced transcript variants have been described.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	AOA
<b>Immunogen:</b>	Fusion protein of human AOA
<b>Full name:</b>	acyloxyacyl hydrolase

<b>SwissProt:</b>	P28039
<b>ELISA Recommended dilution :</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	20-100
<b>WB Predicted band size:</b>	65 kDa
<b>WB Positive control:</b>	HL-60 cell lysate
<b>WB Recommended dilution:</b>	200-1000



