

兔抗 AOC1 多克隆抗体

中文名称：兔抗 AOC1 多克隆抗体

英文名称： Anti-AOC1 rabbit polyclonal antibody

别名： ABP; DAO; KAO; ABP1; DAO1

抗原： AOC1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a metal-binding membrane glycoprotein that oxidatively deaminates putrescine, histamine, and related compounds. The encoded protein is inhibited by amiloride, a diuretic that acts by closing epithelial sodium ion channels. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
Applications:	ELISA, IHC
Name of antibody:	AOC1
Immunogen:	Fusion protein of human AOC1
Full name:	amine oxidase, copper containing 1

Synonyms :	ABP; DAO; KAO; ABP1; DAO1
SwissProt:	P19801
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human colon cancer and Human thyroid cancer
IHC Recommend dilution:	15-50

