

兔抗 AP2A1 多克隆抗体

中文名称: 兔抗 AP2A1 多克隆抗体

英文名称: Anti-AP2A1 rabbit polyclonal antibody

别名: ADTAA; CLAPA1; AP2-ALPHA

抗原: AP2A1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

标记物: Unconjugate

相关类别: 一抗

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes the alpha 1 adaptin subunit of the adaptor protein 2 (AP-2) complex found in clathrin coated vesicles. The AP-2 complex is a heterotetramer consisting of two large adaptins (alpha or beta), a medium adaptin (mu), and a small adaptin (sigma). The complex is part of the protein coat on the cytoplasmic face of coated vesicles which links clathrin to receptors in vesicles. Alternative splicing of this gene results in two transcript variants encoding two different isoforms.
Applications:	ELISA, WB, IHC
Name of antibody:	AP2A1

Immunogen:	Synthetic peptide of human AP2A1
Full name:	adaptor-related protein complex 2, alpha 1 subunit
Synonyms :	ADTAA; CLAPA1; AP2-ALPHA
SwissProt:	O95782
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	108 kDa
WB Positive control:	A431, 231 and NIH/3T3 cell
WB Recommended dilution:	500-2000



