

## 兔抗 AP2B1 多克隆抗体

中文名称：兔抗 AP2B1 多克隆抗体

英文名称： Anti-AP2B1 rabbit polyclonal antibody

别名： adaptor related protein complex 2 beta 1 subunit; ADTB2; AP105B; CLAPB1; AP2-BETA

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： AP2B1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is one of two large chain components of the assembly protein complex 2, which serves to link clathrin to receptors in coated vesicles. The encoded protein is found on the cytoplasmic face of coated vesicles in the plasma membrane. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	AP2B1
<b>Immunogen:</b>	Fusion protein of human AP2B1

<b>Full name:</b>	adaptor related protein complex 2 beta 1 subunit
<b>Synonyms:</b>	ADTB2; AP105B; CLAPB1; AP2-BETA
<b>SwissProt:</b>	P63010
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	105 kDa
<b>WB Positive control:</b>	Mouse brain tissue lysate
<b>WB Recommended dilution:</b>	200-1000

