

兔抗 CLEC4D 多克隆抗体

中文名称：兔抗 CLEC4D 多克隆抗体

英文名称： Anti-CLEC4D rabbit polyclonal antibody

别名： MCL, MPCL, CLEC6, CLEC-6, CLECSF8

抗原： CLEC4D

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. Members of this family share a common protein fold and have diverse functions, such as cell adhesion, cell-cell signaling, glycoprotein turnover, and roles in inflammation and immune response. This gene is closely linked to other CTL/CTLD superfamily members on chromosome 12p13 in the natural killer gene complex region.
Applications:	ELISA, WB
Name of antibody:	CLEC4D
Immunogen:	Fusion protein of human CLEC4D

Full name:	C-type lectin domain family 4, member D
Synonyms :	MCL, MPCL, CLEC6, CLEC-6, CLECSF8
SwissProt:	Q8WXI8
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	25 kDa
WB Positive control:	293T cells
WB Recommended dilution:	500-2000

