

## 兔抗 CLEC6A 多克隆抗体

中文名称：兔抗 CLEC6A 多克隆抗体

英文名称：Anti-CLEC6A rabbit polyclonal antibody

别名：CLEC4N; CLECSF10

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CLEC6A

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Like dectin-1 (CLEC7A; MIM 606264), CLEC6A, or dectin-2, is a type II membrane receptor with an extracellular C-type lectin-like domain fold. However, unlike dectin-1, dectin-2 lacks an immunoreceptor tyrosine-based activation motif (ITAM) in its cytoplasmic domain
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CLEC6A
<b>Immunogen:</b>	Synthetic peptide of human CLEC6A
<b>Full name:</b>	C-type lectin domain family 6, member A
<b>Synonyms :</b>	CLEC4N; CLECSF10
<b>SwissProt:</b>	Q6EIG7

<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

