

兔抗 API5 多克隆抗体

中文名称：兔抗 API5 多克隆抗体

英文名称： Anti-API5 rabbit polyclonal antibody

别名： AAC11; AAC-11

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： API5

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an apoptosis inhibitory protein whose expression prevents apoptosis after growth factor deprivation. This protein suppresses the transcription factor E2F1-induced apoptosis and also interacts with, and negatively regulates Acinus, a nuclear factor involved in apoptotic DNA fragmentation. Its depletion enhances the cytotoxic action of the chemotherapeutic drugs. Multiple alternatively spliced transcript variants encoding different isoforms have been identified.
Applications:	ELISA, WB
Name of antibody:	API5

Immunogen:	Synthetic peptide of human API5
Full name:	apoptosis inhibitor 5
Synonyms :	AAC11; AAC-11
SwissProt:	Q9BZZ5
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	59 kDa
WB Positive control:	293T, A549 and Hela cells
WB Recommended dilution:	500-2000

