

兔抗 APIP 多克隆抗体

中文名称：兔抗 APIP 多克隆抗体

英文名称： Anti-APIP rabbit polyclonal antibody

别名： APIP2, CGI29, CGI-29, MMRP19, dJ179L10.2

抗原： APIP

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	APIP is an APAF1 (MIM 602233)-interacting protein that acts as a negative regulator of ischemic/hypoxic injury. The protein highly expressed in skeletal muscle and heart and binds to the CARD of Apaf-1 in competition with caspase-9. Exogenous APIP inhibits cytochrome c-induced activation of caspase-3 and caspase-9, and suppresses cell death triggered by mitochondrial apoptotic stimuli through inhibiting the downstream activity of cytochrome c released from mitochondria.
Applications:	ELISA, WB, IHC
Name of antibody:	APIP

Immunogen:	Fusion protein of human APIP
Full name:	Probable methylthioribulose-1-phosphate dehydratase
Synonyms :	APIP2, CGI29, CGI-29, MMRP19, dJ179L10.2
SwissProt:	Q96GX9
ELISA Recommended dilution:	500-5000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	27 kDa
WB Positive control:	Mouse heart tissue
WB Recommended dilution:	200-1000



