

兔抗 APOBEC4 多克隆抗体

中文名称：兔抗 APOBEC4 多克隆抗体

英文名称：Anti-APOBEC4 rabbit polyclonal antibody

别名：C1orf169

相关类别：一抗

储存：冷冻（-20℃）

抗原：APOBEC4

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the AID/APOBEC family of polynucleotide (deoxy)cytidine deaminases, which convert cytidine to uridine. Other AID/APOBEC family members are involved in mRNA editing, somatic hypermutation and recombination of immunoglobulin genes, and in innate immunity to retroviral infection.
Applications:	ELISA, IHC
Name of antibody:	APOBEC4
Immunogen:	Fusion protein of human APOBEC4
Full name:	apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 4 (putative)

Synonyms :	C1orf169
SwissProt:	Q8WW27
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

