

兔抗 APOC1 多克隆抗体

中文名称：兔抗 APOC1 多克隆抗体

英文名称：Anti-APOC1 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

抗 原：APOC1

宿 主：Rabbit

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the apolipoprotein C1 family. This gene is expressed primarily in the liver, and it is activated when monocytes differentiate into macrophages. A pseudogene of this gene is located 4 kb downstream in the same orientation, on the same chromosome. This gene is mapped to chromosome 19, where it resides within a apolipoprotein gene cluster. Alternatively spliced transcript variants have been found for this gene, but the biological validity of some variants has not been determined.
Applications:	ELISA, WB
Name of antibody:	APOC1

Immunogen:	Fusion protein of human APOC1
Full name:	apolipoprotein C-I
SwissProt:	P02654
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	9 kDa
WB Positive control:	Human fetal liver tissue
WB Recommended dilution:	500-2000

