

## 兔抗 APOC4 多克隆抗体

中文名称：兔抗 APOC4 多克隆抗体

英文名称： Anti-APOC4 rabbit polyclonal antibody

别名： APO-CIV; APOC-IV

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： APOC4

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a lipid-binding protein belonging to the apolipoprotein gene family. The protein is thought to play a role in lipid metabolism. Polymorphisms in this gene may influence circulating lipid levels and may be associated with coronary artery disease risk. This gene is present in a cluster with other related apolipoprotein genes on chromosome 19. Naturally occurring read-through transcription exists between this gene and the neighboring downstream apolipoprotein C-II (APOC2) gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	APOC4

<b>Immunogen:</b>	Synthetic peptide of human APOC4
<b>Full name:</b>	apolipoprotein C-IV
<b>Synonyms :</b>	APO-CIV; APOC-IV
<b>SwissProt:</b>	P55056
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer and Human brain
<b>IHC Recommend dilution:</b>	15-50

