

兔抗 CLNS1A 多克隆抗体

中文名称：兔抗 CLNS1A 多克隆抗体

英文名称： Anti-CLNS1A rabbit polyclonal antibody

别名： CLCI; ICl_n; CLNS1B

抗原： CLNS1A

储存： 冷冻（-20℃） 避光

宿主： Rabbit

克隆类型： Unconjugate

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

技术规格

Background:	This gene encodes a protein that functions in multiple regulatory pathways. The encoded protein complexes with numerous cytosolic proteins and performs diverse functions including regulation of small nuclear ribonucleoprotein biosynthesis, platelet activation and cytoskeletal organization. The protein is also found associated with the plasma membrane where it functions as a chloride current regulator.
Applications:	WB
Name of antibody:	CLNS1A
Immunogen:	Synthesized peptide derived from C-terminal of human CLNS1A.

Full name:	chloride channel, nucleotide-sensitive, 1A
Synonyms :	CLCI; ICl _n ; CLNS1B
SwissProt:	P54105
WB Predicted band size:	26 kDa
WB Positive control:	Huvec cells, K562 cells, Hela cells, HepG2 cells, A549 cells and Jurkat cells lysates
WB Recommended dilution:	500-3000

