

## 兔抗 APOF 多克隆抗体

中文名称：兔抗 APOF 多克隆抗体

英文名称： Anti-APOF rabbit polyclonal antibody

别名： LTIP; Apo-F

抗原： APOF

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	The product of this gene is one of the minor apolipoproteins found in plasma. This protein forms complexes with lipoproteins and may be involved in transport and/or esterification of cholesterol.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	APOF
<b>Immunogen:</b>	Synthesized peptide derived from internal of human APOF.
<b>Full name:</b>	apolipoprotein F
<b>Synonyms :</b>	LTIP; Apo-F
<b>SwissProt:</b>	Q13790

<b>WB Predicted band size:</b>	35 kDa
<b>WB Positive control:</b>	A549 cells, COLO cells, Jurkat cells and Hela cells lysates
<b>WB Recommended dilution:</b>	500-3000

