

## 兔抗 APOL1 多克隆抗体

中文名称：兔抗 APOL1 多克隆抗体

英文名称： Anti-APOL1 rabbit polyclonal antibody

别名： APOL; APO-L; FSGS4; APOL-I

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： APOL1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

This gene encodes a secreted high density lipoprotein which binds to apolipoprotein A-I. Apolipoprotein A-I is a relatively abundant plasma protein and is the major apoprotein of HDL. It is involved in the formation of most cholesteryl esters in plasma and also promotes efflux of cholesterol from cells. This apolipoprotein L family member may play a role in lipid exchange and transport throughout the body, as well as in reverse cholesterol transport from peripheral cells to the liver. Several different transcript variants encoding different isoforms have been found for this gene.

<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	APOL1
<b>Immunogen:</b>	Synthetic peptide of human APOL1
<b>Full name:</b>	apolipoprotein L, 1
<b>Synonyms:</b>	APOL; APO-L; FSGS4; APOL-I
<b>SwissProt:</b>	O14791
<b>IHC positive control:</b>	Human gastric cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	44 kDa
<b>WB Positive control:</b>	Human plasma tissue
<b>WB Recommended dilution:</b>	200-1000



