

兔抗 CLUAP1 多克隆抗体

中文名称：兔抗 CLUAP1 多克隆抗体

英文名称： Anti-CLUAP1 rabbit polyclonal antibody

别名： FAP22

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： CLUAP1

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene contains a single coiled-coil region. Alternative splicing results in multiple transcript variants and protein isoforms.
Applications:	ELISA, IHC
Name of antibody:	CLUAP1
Immunogen:	Fusion protein of human CLUAP1
Full name:	clusterin associated protein 1
Synonyms :	FAP22
SwissProt:	Q96AJ1
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human cervical cancer and human tonsil

IHC Recommend dilution:

100-300

