

兔抗 ARCNI 多克隆抗体

中文名称：兔抗 ARCNI 多克隆抗体

英文名称： Anti-ARCNI rabbit polyclonal antibody

别名： archain 1; COPD

相关类别： 一抗

抗原： ARCNI

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene maps in a region, which include the mixed lineage leukemia and Friend leukemia virus integration 1 genes, where multiple disease-associated chromosome translocations occur. It is an intracellular protein. Archain sequences are well conserved among eukaryotes and this protein may play a fundamental role in eukaryotic cell biology. It has similarities to heat shock proteins and clathrin-associated proteins, and may be involved in vesicle structure or trafficking.
Applications:	ELISA, IHC
Name of antibody:	ARCNI

Immunogen:	Synthetic peptide of human ARCN1
Full name:	archain 1
Synonyms:	COPD
SwissProt:	P48444
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

