

兔抗 CLU/APO J 多克隆抗体

中文名称：兔抗 CLU/APO J 多克隆抗体

英文名称：Anti-CLU/APO J rabbit polyclonal antibody

别名：CLI; AAG4; APOJ; CLU1; CLU2; KUB1; SGP2; APO-J; SGP-2; SP-40; TRPM2; TRPM-2; NA1/NA2

抗原：CLU/APO J

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Rabbit

反应种属：Human, Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	The protein encoded by this gene is a secreted chaperone that can under some stress conditions also be found in the cell cytosol. It has been suggested to be involved in several basic biological events such as cell death, tumor progression, and neurodegenerative disorders. Alternative splicing results in both coding and non-coding variants.
Applications:	ELISA, IHC
Name of antibody:	CLU/APO J
Immunogen:	Fusion protein of human CLU/APO J

Full name:	clusterin
Synonyms :	CLI; AAG4; APOJ; CLU1; CLU2; KUB1; SGP2; APO-J; SGP-2 ; SP-40; TRPM2; TRPM-2; NA1/NA2
SwissProt:	P10909
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain and human liver cancer
IHC Recommend dilution:	50-350



