

## 兔抗 ARF6 多克隆抗体

- 中文名称: 兔抗 ARF6 多克隆抗体
- 英文名称: Anti-ARF6 rabbit polyclonal antibody
- 储存: 冷冻(-20℃)
- 抗原: ARF6
- 宿 主: Rabbit
- 相关类别: 一抗
- 反应种属: Human, Mouse, Rat
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

## 技术规格

Background:	This gene encodes a member of the human ARF gene family, which is part of the RAS superfamily. The ARF g enes encode small guanine nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of chol era toxin and play a role in vesicular trafficking and as activators of phospholipase D. The product of this gene is localized to the plasma membrane, and regulates ves icular trafficking, remodelling of membrane lipids, and s ignaling pathways that lead to actin remodeling. A pse udogene of this gene is located on chromosome 7.
Applications:	ELISA, WB, IHC
Name of antibody:	ARF6
Immunogen:	Synthetic peptide of human ARF6



Full name:	ADP-ribosylation factor 6
SwissProt:	P62330
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human colon cancer
IHC Recommend dilution:	10-50
WB Predicted band size:	20 kDa
WB Positive control:	HepG2 and 231 cells
WB Recommended dilution:	500-2000





