

兔抗 ARHGAP35(Phospho-Tyr1105) 多克隆 抗体

中文名称：兔抗 ARHGAP35(Phospho-Tyr1105) 多克隆抗体

英文名称： Anti-ARHGAP35(Phospho-Tyr1105) rabbit polyclonal antibody

别名： GRF-1; GRLF1; P190A; P190-A; p190RhoGAP; p190ARhoGAP

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： ARHGAP35(Phospho-Tyr1105)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The human glucocorticoid receptor DNA binding factor, which associates with the promoter region of the glucocorticoid receptor gene (hGR gene), is a repressor of glucocorticoid receptor transcription. The amino acid sequence deduced from the cDNA sequences show the presence of three sequence motifs characteristic of a zinc finger and one motif suggestive of a leucine zipper in which 1 cysteine is found instead of all leucines.
Applications:	WB

Name of antibody:	ARHGAP35(Phospho-Tyr1105)
Immunogen:	Peptide sequence around phosphorylation site of Tyrosine 1105 (N-I-Y(p)-S-V) derived from Human GRF-1.
Full name:	Rho GTPase activating protein 35
Synonyms :	GRF-1; GRLF1; P190A; P190-A; p190RhoGAP; p190ARhoGAP
SwissProt:	Q9NRY4
WB Predicted band size:	171 kDa
WB Positive control:	Mouse brain tissue lysates
WB Recommended dilution:	500-1000

