

## 兔抗 ARHGAP45 多克隆抗体

- 中文名称：兔抗 ARHGAP45 多克隆抗体
- 英文名称：Anti-ARHGAP45 rabbit polyclonal antibody
- 别名：Rho GTPase activating protein 45; HA-1; HMHA1; HLA-HA1
- 抗原：ARHGAP45
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Contains a GTPase activator for the Rho-type GTPases (RhoGAP) domain that would be able to negatively regulate the actin cytoskeleton as well as cell spreading. However, also contains N-terminally a BAR-domain which is able to play an autoinhibitory effect on this RhoGAP activity.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	ARHGAP45
<b>Immunogen:</b>	Synthetic peptide of human ARHGAP45
<b>Full name:</b>	Rho GTPase activating protein 45
<b>Synonyms:</b>	HA-1; HMHA1; HLA-HA1
<b>SwissProt:</b>	Q92619

<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	125 kDa
<b>WB Positive control:</b>	A375 cell lysate
<b>WB Recommended dilution:</b>	500-2000

