

兔抗 CMTM4 多克隆抗体

中文名称：兔抗 CMTM4 多克隆抗体

英文名称：Anti-CMTM4 rabbit polyclonal antibody

别名：CKLF like MARVEL transmembrane domain containing 4; CKLFSF4

相关类别：一抗

抗原：CMTM4

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and the transmembrane 4 superfamilies of signaling molecules. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 16. Alternatively spliced transcript variants encoding different isoforms have been identified.
Applications:	ELISA, IHC
Name of antibody:	CMTM4
Immunogen:	Synthetic peptide of human CMTM4
Full name:	CKLF like MARVEL transmembrane domain containing 4

Synonyms:	CKLFSF4
SwissProt:	Q8IZR5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human tonsil
IHC Recommend dilution:	50-200



